

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-HSP40/DNAJB1 Polyclonal antibody
PRODUCT REFERENCE	ARO-A14044
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human HSP40 / DNAJB1 (Met1-Ile340).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-DnaJ protein homolog 1,HSP40,DnaJ homolog subfamily B member 1,DNAJB1,DNAJ1,Heat shock 40 kDa protein 1,hDj-1,HDJ1,Human DnaJ protein 1,Heat shock protein 40,HSPF1,