

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|--|
| PRODUCT NAME | Anti-NFKB1 Polyclonal antibody |
| PRODUCT REFERENCE | ARO-A14076 |
| CLONALITY | Polyclonal Antibody |
| ISOTYPE | IgG |
| TARGET | E. coli - derived recombinant Human NFKB1 (Met1-Ile969). |

PRODUCT INFORMATION

| | |
|--------------------------------|---|
| ORIGIN SPECIES | <i>Oryctolagus cuniculus (Rabbit)</i> |
| FORM | Liquid |
| BUFFER | 0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300 |
| STABILITY & STORAGE | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| ALIAS | Anti-Nuclear factor NF-kappa-B p105 subunit,EBP-1,NFKB1,Nuclear factor of kappa light polypeptide gene enhancer in B-cells 1,DNA-binding factor KBF1, |