

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-POU2F1 Polyclonal antibody
PRODUCT REFERENCE	ARO-A14112
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human POU2F1 (Ala413-Val671).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Oct-1,OTF-1,NF-A1,Octamer-binding transcription factor 1,Octamer-binding protein 1,POU2F1,OCT1,POU domain,class 2,transcription factor 1,OTF1,