

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-CD331/FGFR1 Polyclonal antibody
PRODUCT REFERENCE	ARO-A14140
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human CD331 / FGFR1 (Trp471-Arg821).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-FLT2,HBGFR,BFGFR,FGFR-1,CEK,Fibroblast growth factor receptor 1,FGFR1,bFGF-R-1,Fms-like tyrosine kinase 2,FGFBR,FLT-2,CD331,Basic fibroblast growth factor receptor 1,FLG,Proto-oncogene c-Fgr,N-sam,