

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-CXCL1/NAP-3/GRO-alpha Polyclonal antibody
PRODUCT REFERENCE	ARO-A14164
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human CXCL1 / NAP-3 / GRO-alpha (Ala35-Asn107).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-CXCL1,GROA,GRO1,Neutrophil-activating protein 3,GRO,C-X-C motif chemokine 1,NAP-3,Growth-regulated alpha protein,MGSA,SCYB1,Melanoma growth stimulatory activity,GRO-alpha(1-73),