

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-FGF2 Polyclonal antibody
PRODUCT REFERENCE	ARO-A14173
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human FGF2 (Pro143-Ser288(Cys211Ser, Cys229Ser)).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-bFGF,FGF2,Basic fibroblast growth factor,HBGF-2,FGFB,FGF-2,Heparin-binding growth factor 2,Fibroblast growth factor 2,