

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-HEXB Polyclonal antibody
PRODUCT REFERENCE	ARO-A14190
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human HEXB (Lys315-Tyr456).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-N-Acetyl-Î²-glucosaminidase (NAG),chitobiase,NAG,Beta-hexosaminidase subunit beta,Beta-N-acetylhexosaminidase subunit beta,N-acetyl-beta-glucosaminidase subunit beta,Cervical cancer proto-oncogene 7 protein,Hexosaminidase subunit B,HCC-7,HEXB,