

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-ENO1 Polyclonal antibody
PRODUCT REFERENCE	ARO-A14207
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human ENO1 (Pro267-Leu354).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Non-neural enolase,2-phospho-D-glycerate hydro-lyase,MPB-1,ENO1,Plasminogen-binding protein,Enolase 1,C-myc promoter-binding protein,Alpha-enolase,MBP-1,MBPB1,MPB1,ENO1L1,NNE,Phosphopyruvate hydratase,