

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-KNG1 Polyclonal antibody
PRODUCT REFERENCE	ARO-A14259
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human KNG1 (Arg58-Leu176).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Kininogen-1, Kallidin II, Fitzgerald factor, KNG, Alpha-2-thiol proteinase inhibitor, Kallidin I, Williams-Fitzgerald-Flaujeac factor, High molecular weight kininogen, HMWK, KNG1, Ile-Ser-Bradykinin, BDK,