

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-PRDM1 Polyclonal antibody
PRODUCT REFERENCE	ARO-A14311
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human PRDM1 (His493-Lys678).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-BLIMP-1,PRDI-BF1,PRDI-binding factor 1,BLIMP1,Positive regulatory domain I-binding factor 1,Beta-interferon gene positive regulatory domain I-binding factor,PR domain-containing protein 1,PR domain zinc finger protein 1,PRDM1,