

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-MACROH2A1 Polyclonal antibody
PRODUCT REFERENCE	ARO-A14317
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human MACROH2A1 (Ile109-Leu369).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Histone macroH2A1,mH2A1,H2A/y,MACROH2A1,Medulloblastoma antigen MU-MB-50.205,Core histone macro-H2A.1,H2AFY,Histone H2A.y,