

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-CCL16/HCC-4 Polyclonal antibody
<b>PRODUCT REFERENCE</b>	ARO-A14339
<b>CLONALITY</b>	Polyclonal Antibody
<b>ISOTYPE</b>	IgG
<b>TARGET</b>	E. coli - derived recombinant Human CCL16 / HCC-4 (Gln24-Gln120).

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Oryctolagus cuniculus (Rabbit)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	Anti-Chemokine CC-4,IL-10-inducible chemokine,NCC-4,SCYA16,Liver-expressed chemokine,MTN-1,LMC,Monotactin-1,NCC4,Small-inducible cytokine A16,Chemokine LEC,CCL16,C-C motif chemokine 16,ILINCK,HCC-4,Lymphocyte and monocyte chemoattractant,LCC-1,