

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-P4HA2 Polyclonal antibody
PRODUCT REFERENCE	ARO-A14341
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human P4HA2 (Ser207-Glu520).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-P4HA2,4-PH alpha-2,Procollagen-proline,2-oxoglutarate-4-dioxygenase subunit alpha-2,Prolyl 4-hydroxylase subunit alpha-2,