

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-CCN1 Polyclonal antibody
PRODUCT REFERENCE	ARO-A14361
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human CCN1 (Leu178-Asp381).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Insulin-like growth factor-binding protein 10,Protein GIG1,IGF-binding protein 10,CYR61,Cellular communication network factor 1,IBP-10,Cysteine-rich angiogenic inducer 61,GIG1,Protein CYR61,CCN1,IGFBP-10,IGFBP10,CCN family member 1,