

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-DRP1/DNM1L Polyclonal antibody
PRODUCT REFERENCE	ARO-A14365
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human DRP1 / DNM1L (Met1-Arg516).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Dynamin family member proline-rich carboxyl-terminal domain less,DVLP,Dynamin-like protein IV,Dynamin-related protein 1,Dynamin-like protein,DLP1,Dymple,Dynamin-1-like protein,HdynIV,Dnm1p/Vps1p-like protein,DRP1,DNM1L,Dynamin-like protein 4,