

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-IGF2BP3 Polyclonal antibody
PRODUCT REFERENCE	ARO-A14367
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human IGF2BP3 (Asn2-Ile343).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Insulin-like growth factor 2 mRNA-binding protein 3,VICKZ3,IGF2BP3,IMP3,KH domain-containing protein overexpressed in cancer,IMP-3,VICKZ family member 3,IGF-II mRNA-binding protein 3,KOC1,hKOC,IGF2 mRNA-binding protein 3,