

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Ganglioside GD3 Antibody (R24)
PRODUCT REFERENCE	ARO-A14592
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	Disialoganglioside GD3, GD3 ganglioside, $\hat{1}\pm$ -N-acetylneuraminosyl-(2 $\hat{2}$ '8)- $\hat{1}\pm$ -N-acetylneuraminosyl-(2 $\hat{2}$ '3)- $\hat{1}^2$ -D-galactosyl-(1 $\hat{1}$ '4)- $\hat{1}^2$ -D-glucosyl-(1 $\hat{1}$ '1')-ceramide

PRODUCT INFORMATION

ORIGIN SPECIES	<i>General</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS