



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Staphylococcus aureus isdB Antibody (SAA1456)
PRODUCT REFERENCE	ARO-A14597
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a
TARGET	isdB, frpB, sasJ, sirH, Iron-regulated surface determinant protein B, Fur-regulated protein B, Staphylococcal iron-regulated protein H, Staphylococcus aureus surface protein J

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Staphylococcus aureus</i> (strain MSSA476) (Q6GA86), <i>Staphylococcus aureus</i> (strain bovine RF122 / ET3-1) (Q2YX96), <i>Staphylococcus aureus</i> (strain COL) (Q5HGV5)
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY