

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Staphylococcus aureus isdA/isdB Antibody (SAA1457)
PRODUCT REFERENCE	ARO-A14598
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a
TARGET	isdA, frpA, stbA, Iron-regulated surface determinant protein A, Fur-regulated protein A, Staphylococcal transferrin-binding protein A, isdB, frpB, sasJ, sirH, Iron-regulated surface determinant protein B, Fur-regulated protein B, Staphylococcal iron-regul

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Staphylococcus aureus</i> (strain NCTC 8325 / PS 47) (Q2FZE9), <i>Staphylococcus aureus</i> (strain NCTC 8325 / PS 47) (Q2FZF0)
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	