

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Bacillus anthracis PA/pagA/PA-83 Antibody (SAA1376)
PRODUCT REFERENCE	ARO-A14617
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Protective antigen, PA, Anthrax toxins translocating protein, PA-83, PA83, Protective antigen PA-20, PA-20, PA20, Protective antigen PA-63, PA-63, PA63, pagA, pag, pXO1-110, BXA0164, GBAA_pXO1_0164

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Bacillus anthracis</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS