



ProteoGenix

ProteoGenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

ProteoGenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human KRT8/18 Antibody (SAA1385)
PRODUCT REFERENCE	ARO-A14633
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Type-II keratin Kb8, CYK8, Keratin-8, KRT8, Keratin, type II cytoskeletal 8, Cytokeratin-8, CK-8, K8, Cytokeratin-18, CYK18, K18, Keratin, type I cytoskeletal 18, KRT18, CK-18, Keratin-18, Cell proliferation-inducing gene 46 protein

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY