



ProteoGenix

ProteoGenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

ProteoGenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Gluten Antibody (SAA0755)
PRODUCT REFERENCE	ARO-A14751
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	Gluten (prolamin and glutenin) from wheat (<i>Triticum aestivum</i>), barley (<i>Hordeum vulgare</i>) and rye (<i>Secale cereale</i>)

PRODUCT INFORMATION

ORIGIN SPECIES	<i>wheat (Triticum aestivum)</i> , <i>barley (Hordeum vulgare)</i> , <i>rye (Secale cereale)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY