

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-mitF/aspF1 Antibody (SAA0653)
PRODUCT REFERENCE	ARO-A14784
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	mitF, aspF1, Ribonuclease mitogillin, Allergen Asp f I, Allergen I/a, IgE-binding ribotoxin, Major allergen Asp f 1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Neosartorya fumigata</i> (strain ATCC MYA-4609 / Af293 / CBS 101355 / FGSCA1100) (<i>Aspergillus fumigatus</i>)
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS