

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Bet v 1/Allergen Bet v I-A Antibody (SAA0669)
PRODUCT REFERENCE	ARO-A14785
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	BETVIA, BETVI, Major pollen allergen Bet v 1-A, Allergen Bet v 1-A, ALNGI, Major pollen allergen Aln g 1, Allergen Aln g I, CORAI, Major pollen allergen Cor a 1, Allergen Cor a I, MALDI, MALD1, Major allergen Mal d 1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Betula pendula</i> (European white birch) (<i>Betula verrucosa</i>)
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS