

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-PHLPVa/Allergen Phl p Va Antibody (SAA0702)
PRODUCT REFERENCE	ARO-A14791
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	PHLPVa, Pollen allergen Phl p 5a, Allergen Phl p Va

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Phleum pratense (Common timothy)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS