

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Mal f 1/Allergen Pit o 1 Antibody (SAA0647)
PRODUCT REFERENCE	ARO-A14821
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-hFc
TARGET	Major allergen Mal f 1, Allergen Pit o 1, Mal f 1, Cell surface protein Malf1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Malassezia furfur (Pityriasis versicolor infection agent) (Pityrosporumfurfur)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS