



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-LPS/Lipopolysaccharide Antibody (SAA0585)
PRODUCT REFERENCE	ARO-A14880
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a
TARGET	Lipopolysaccharide (LPS), O-polysaccharide moiety (Klebsiella pneumoniae O1 serotype O-antigen)Galactan-11, D-galactan-II, D-gal II, gal-11, [α -3)- β -d-Galp-(1 α -3)- β -d-Galp-(1 α -3)

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Klebsiella pneumoniae</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY