



# ProteoGenix

**ProteoGenix S.A.S**

19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

[products@proteogenix.fr](mailto:products@proteogenix.fr)

**ProteoGenix INC.**

1000 Bearcat Way, Suite  
100, Morrisville, NC 27560,  
USA

Tel +1(919)234-1277

[products@proteogenix.fr](mailto:products@proteogenix.fr)

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-LPS/Lipopolysaccharide Antibody (SAA0577)
<b>PRODUCT REFERENCE</b>	ARO-A14907
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG2a
<b>TARGET</b>	Lipopolysaccharide (LPS), O-polysaccharide moiety ( <i>Pseudomonas aeruginosa</i> serogroup O6 O-polysaccharide)

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Pseudomonas aeruginosa</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

### ALIAS

You can contact us at [products@proteogenix.fr](mailto:products@proteogenix.fr) for any question

**OUR PRODUCTS ARE FOR RESEARCH USE ONLY**