

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human ADAM15/MDC15 Antibody (SAA0494)
PRODUCT REFERENCE	ARO-A14940
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Disintegrin and metalloproteinase domain-containing protein 15, Metalloproteinase-like, disintegrin-like, and cysteine-rich protein 15, ADAM15, MDC15, Metalloprotease RGD disintegrin protein, Metargidin, MDC-15, ADAM 15

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS