

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Vaccinia virus/VACV A27L Antibody (1G6)
PRODUCT REFERENCE	ARO-A14945
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Protein A27, VACWR150, A27L, Vaccinia virus (VACV)

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Vaccinia virus (strain Western Reserve) (VACV) (Vaccinia virus (strainWR))</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS