



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-FMDV Genome polyprotein Antibody (201#)
PRODUCT REFERENCE	ARO-A14961
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Genome polyprotein, Leader protease, Lpro, 3.4.22.46, Capsid protein VP0, VP4-VP2, Capsid protein VP4, P1A, Virion protein 4, Capsid protein VP2, P1B, Virion protein 2, Capsid protein VP3, P1C, Virion protein 3, Capsid protein VP1, P1D, Virion protein 1,

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Foot-and-mouth disease virus - type O, Foot-and-mouth disease virus (isolate Bovine/Brazil/A24Cruzeiro/1955serotype A) (FMDV)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY