

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-EBV/HHV-4 LMP1/BNLF1 Antibody (B8)
PRODUCT REFERENCE	ARO-A14971
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Latent membrane protein 1, LMP-1, Protein p63, Protein p25, LMP1, BNLF1, Epstein-Barr virus (EBV), Human herpesvirus 4 (HHV-4), HHV4

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Epstein-Barr virus (strain B95-8) (HHV-4) (Human herpesvirus 4)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS