

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-EBV/HHV-4 LMP1/BNLF1 Antibody (A3H5)
<b>PRODUCT REFERENCE</b>	ARO-A15009
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1, kappa
<b>TARGET</b>	Latent membrane protein 1, LMP-1, Protein p63, Protein p25, LMP1, BNLF1, Epstein-Barr virus (EBV), Human herpesvirus 4 (HHV-4), HHV4

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Epstein-Barr virus (strain B95-8) (HHV-4) (Human herpesvirus 4)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

### ALIAS