



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-EBV/HHV-4 EBNA1/BKRF1 Antibody (3D4)
PRODUCT REFERENCE	ARO-A15010
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Epstein-Barr nuclear antigen 1, EBNA-1, EBV nuclear antigen 1, EBNA1, BKRF1, Epstein-Barr virus (EBV), Human herpesvirus 4 (HHV-4), HHV4

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Epstein-Barr virus (strain B95-8) (HHV-4) (Human herpesvirus 4)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY