

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-Phospho-CREB1 (pS119) Antibody (P5.G10)
<b>PRODUCT REFERENCE</b>	ARO-A15039
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG2a, kappa
<b>TARGET</b>	cAMP-responsive element-binding protein 1, CREB-1, CREB1, Cyclic AMP-responsive element-binding protein 1

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Human, Mouse, Rat, Bovine</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

### ALIAS