

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Phospho- $\hat{\text{I}}^{\text{s}}$ -H2AX (pS139) Antibody (C6)
PRODUCT REFERENCE	ARO-A15041
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-mFc (IgG2a)
TARGET	H2a/x, Histone H2A.X, H2AFX, Histone H2AX, H2AX, gamma H2A.X (phospho S139)

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS