

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Phospho-POLR2A-CTD (p(YSP TSPSYSP TSPS)-2ndS2) Antibody (2ndS2)
PRODUCT REFERENCE	ARO-A15063
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	DNA-directed RNA polymerase II subunit A, RNA-directed RNA polymerase II subunit RPB1, RNA polymerase II subunit B1, POLR2A, POLR2, DNA-directed RNA polymerase II subunit RPB1, DNA-directed RNA polymerase III largest subunit

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS