



ProteoGenix

ProteoGenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

ProteoGenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-HCV NS1/gp68/gp70/Envelope glycoprotein E2 Antibody (19B3)
PRODUCT REFERENCE	ARO-A15096
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	NS1, gp68, gp70, Envelope glycoprotein E2, Genome polyprotein, Hepatitis C Virus (HCV)

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Hepatitis C virus genotype 1a (isolate H77) (HCV), Hepatitis C virus genotype 1b (isolate Con1) (HCV), Hepatitis C virus genotype 2a (isolate HC-J6) (HCV), Hepatitis C virus S52/JFH1, Hepatitis C virus genotype 4a (isolate ED43) (HCV), Hepatitis C virus g</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY