

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human SMYD2 IP Antibody (SAA0229)
PRODUCT REFERENCE	ARO-A15141
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a
TARGET	HSKM-B, N-lysine methyltransferase SMYD2, Histone methyltransferase SMYD2, SMYD2, SET and MYND domain-containing protein 2, Lysine N-methyltransferase 3C, KMT3C

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS