

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-Human TAF1 IP Antibody (SAA0232)
<b>PRODUCT REFERENCE</b>	ARO-A15147
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG2a
<b>TARGET</b>	p250, TAF(II)250, CCGS, BA2R, TAF1, Transcription initiation factor TFIID subunit 1, Cell cycle gene 1 protein, TAFII-250, TAF2A, CCG1, TBP-associated factor 250 kDa, Transcription initiation factor TFIID 250 kDa subunit, TAFII250

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Homo sapiens (Human)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

#### ALIAS