



# Proteogenix

**Proteogenix S.A.S**

19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

**Proteogenix INC.**

1000 Bearcat Way, Suite  
100, Morrisville, NC 27560,  
USA

Tel +1(919)234-1277

products@proteogenix.fr

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-HCV NS1/gp68/gp70/Envelope glycoprotein E2 Antibody (AP33)
<b>PRODUCT REFERENCE</b>	ARO-A15269
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1, kappa
<b>TARGET</b>	NS1, gp68, gp70, Envelope glycoprotein E2, Genome polyprotein, Hepatitis C Virus (HCV)

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Hepatitis C virus genotype 1a (isolate H77) (HCV) &amp; Hepatitis C virus genotype 1b (isolate BK) (HCV) &amp; Hepatitis C virus (isolate Glasgow) (HCV)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

### ALIAS

You can contact us at [products@proteogenix.fr](mailto:products@proteogenix.fr) for any question

**OUR PRODUCTS ARE FOR RESEARCH USE ONLY**