

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-HCV NS1/gp68/gp70/Envelope glycoprotein E2 Antibody (DAO5)
PRODUCT REFERENCE	ARO-A15288
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	NS1, gp68, gp70, Envelope glycoprotein E2, Genome polyprotein, Hepatitis C Virus (HCV)

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Hepatitis C virus genotype 1b (strain HC-J4) (HCV) & Hepatitis C virus genotype 2a (isolate JFH-1) (HCV)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS