

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-HCV Genome polyprotein Antibody (cVH-E2)
PRODUCT REFERENCE	ARO-A15291
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Genome polyprotein, Hepatitis C Virus (HCV)

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Hepatitis C virus genotype 1b (isolate BK) (HCV)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS