

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-HeV/NiV F/Fusion glycoprotein F0 Antibody (mAb66)
PRODUCT REFERENCE	ARO-A15302
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	F, F ₀ , Fusion glycoprotein F ₀ , Protein F, Nipah Virus (NiV) / Hendra Virus (HeV)

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Nipah virus</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS