

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Recombinant Anti-HLA-A2-peptide (FMNKFIYEI) Complex (AMC-48)
PRODUCT REFERENCE	ARO-A15306
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	HLAA,Human leukocyte antigen A,HLA class I histocompatibility antigen,A alpha chain,HLA-A,Alpha-fetoprotein,Alpha-1-fetoprotein,Alpha-fetoglobulin,AFP,HPAFP

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS