

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-HeV/NiV Glycoprotein G Antibody (HENV-26)
PRODUCT REFERENCE	ARO-A15322
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	G, Glycoprotein G, Attachment glycoprotein, Receptor binding protein, RBP, Nipah Virus (NiV) / Hendra Virus (HeV)

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Hendra virus (isolate Horse/Australia/Hendra/1994) & Nipah virus</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS