



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|---|
| PRODUCT NAME | Recombinant Anti-HLA-A2-peptide (YLEPGPVT) Complex (2B2) |
| PRODUCT REFERENCE | ARO-A15324 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG1 |
| TARGET | HLAA,Human leukocyte antigen A,HLA class I histocompatibility antigen,A alpha chain,HLA-A,ME20M,Melanocyte protein PMEL,Silver locus protein homolog,Secreted melanoma-associated ME20 antigen,D12S53E,95 kDa melanocyte-specific secreted glycoprotein,SILV,P1 |

PRODUCT INFORMATION

| | |
|--------------------------------|--|
| ORIGIN SPECIES | <i>Homo sapiens (Human)</i> |
| FORM | Liquid |
| BUFFER | 0.01M PBS (pH 7.4) |
| STABILITY & STORAGE | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |

ALIAS

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY