

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-HRSV-A F/Fusion glycoprotein F0 Antibody (HNK20-hu10)
PRODUCT REFERENCE	ARO-A15335
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	F, Fusion glycoprotein F0, Protein F, Human Respiratory Syncytial Virus (HRSV),RSV

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Human respiratory syncytial virus A (strain A2)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS