

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human IL12B/IL-12 p40/NKSF2 Antibody (SAA0381)
PRODUCT REFERENCE	ARO-A15574
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	IL12B/IL-12 p40/NKSF2

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 0.09% Sodium azide
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS