

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD85d/LILRB2 Monoclonal Antibody
PRODUCT REFERENCE	ARO-A15773
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	CD85d/LILRB2

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	CD85 antigen-like family member D, Monocyte/macrophage immunoglobulin-like receptor 10, MIR-10, ILT4, LIR-2, ILT-4, Immunoglobulin-like transcript 4, MIR10, CD85d, Leukocyte immunoglobulin-like receptor subfamily B member 2, Leukocyte immunoglobulin-like